

# GET YOUR LEVEL UP WITH US!

Put the rest on mute!



E  
Q  
F  
Y



BY  
earKART

## High Quality Sound Processing Through 16 Channels & 12 Bands



### The Best Your Money can Buy!



Make call &  
stream music



IP65 Water &  
Dust Resistance



Rechargeable  
Hearing Aid

Awesome  
Price

Compare  
Features

Compare  
Performance

eQFY

BY  
earKART

# IP65

## Water & Dust Resistance



Enjoy Rain | Enjoy Sweat | Enjoy Life



BY  
earKART

# Fully customizable for **Personalized sound via app**

(Even a 10 year old can do it, so can you)



Convenience



Adjust sound through app  
yourself anywhere



No frequent  
clinic visit

Peace of mind like never before



# Make Call & stream music at your fingertips

Stream calls and music to both ears with  
**classic bluetooth connectivity** through EQFY

Freedom | Convenience | Comfort





eQFY

BY  
earKART



# 3-in-1

compact power bank



Charges 3 times  
(Power Bank)



UV Disinfection



Portable Case



Light weight | Dependable | Easy

eQFY

BY  
earKART

# Say goodbye to complicated battery changes with EQFY



COMFORT



Battery lasting upto 26 hours.  
The carrying case can also charge  
the device up to 3 times.

# FAQ'S

## **Q1. What is EQFY?**

EQFY is a rechargeable Receiver in the canal (RIC) hearing aid.

## **Q2. Is EQFY an analogue or digital hearing aid?**

EQFY is a digital hearing aid that can be programmed according to the users comfort.

## **Q3. Can I swim while using an EQFY device?**

Though EQFY is water resistant, means it can withstand some water, it is recommended not to use hearing aids while swimming.

## **Q4. Which smartphones are compatible with EQFY hearing aids?**

All smartphones having Bluetooth connectivity option are compatible with EQFY.

## **Q5. How often should I charge my rechargeable hearing aids?**

You should charge your hearing aids every night and there is no need to fully discharge the battery before placing your hearing aids back in the charger

## **Q6. How long will the battery last on a full charge?**

The daily performance of the battery depends on a number of factors such as your hearing loss, the sound environment around you, the age of the battery, and whether you are streaming sound from your phone, TV, or other devices. Below is battery consumption time in different modes

- Hearing Aid Mode: > 20 Hours.
- Bluetooth Mode: > 15 Hours
- Bluetooth Streaming: > 4 Hours
- Bluetooth on Calling: > 2 Hours

## **Q7. How long does it take to fully charge the hearing aids?**

If the battery is fully drained, it takes approximately two hours to fully charge your hearing aids. Quick 15 mins charge provide 8 hours of back-up

## **Q8. How long can I leave my hearing aids in the charger?**

The charger is the best place to store your hearing aids. Any time your hearing aids are not in use, place them in the charger for charging. This will ensure your hearing aids are always charged.

## **Q9. Will over charging damage my phone?**

Your hearing aid will automatically go into power off mode once fully charged.

## **Q10. How do I turn my rechargeable hearing aids on or off ?**

Your hearing aids automatically turn on when you remove them from the charger. The hearing aids automatically turn off when you place them in the charger.

## **Q11. What is the warranty on EQFY Hearing Aids?**

earKart offers two years warranty on EQFY where as one year warranty on portable charger case.



**Q12. Will EQFY help me in streaming music and answering Phone calls?**  
Yes, once connected to Bluetooth & Earkart Hearing App. You can stream music and your phone calls.

**Q13. What type of battery do the rechargeable hearing aids has?**  
Lithium-ion batteries are used in rechargeable hearing aids.

**DOWNLOAD THE "earKART HEARING " APP**



\*The "earKART Hearing" app is compatible with Android and iOS. You can personalize your hearing experience.

Step 1: Download the "earKART Hearing" from the App Store or Google Play Store is free of charge.

Step 2: Enable the Bluetooth on your phone, search the device name of "EQFY"

Step 3: Tap to connect your hearing aid with your phone, you can adjust the volume and customized adjustment.

\*Search and connect the device name of "EQFY-L" or "EQFY-R".



## Experience the best service

earKART is a group of professionals and industry stalwarts who work towards increasing awareness about hearing loss & making the hearing aid selection process transparent, consistent and fair every time.

## Most unique earKART's benefits



**Free Insurance**  
for your Hearing Aid



**50% Instant**  
Wallet Credit



**Exchange/Upgrade**  
Hearing Aid



**Strong**  
Customer Support



**Free Senior Audiologist**  
Consultation

**Exclusive benefits Upto**

**₹25,000**



# Most Economical Rechargable hearing aid with Top featues.



earKart hearing app at your fingertips.



Binaural Bluetooth Connectivity



Rechargeable Hearing aid



IP65 Water & Dust Resistance



16 Channels and 12 Bands



Discreet Look & Design



Adjust sound through app yourself



3-in-1 compact carrying case



Self assessment



Noise Reduction



+91 9289097578



info@eartkart.in



www.eartkart.com



Eartkart Pvt. Ltd.  
A-133, Block A, Sec 63,  
Noida-201301, India